

NGS data formats and Quality Control

Jacqui Keane

@drjkeane

drjkeane@gmail.com

Adapted from slides provided by Petr Danecek
petr.danecek@sanger.ac.uk

FASTQ

- | Unaligned read sequences with base qualities

SAM/BAM


- | Unaligned or aligned reads
- | Text and binary formats

CRAM

- | Better compression than BAM

VCF/BCF

- | Flexible variant call format
- | Arbitrary types of sequence variation
- | SNPs, indels, structural variations



Specifications maintained by the Global Alliance for Genomics and Health

FASTA - reference genome

```
>1 dna:chromosome chromosome:GRCh37:1:1:249250621:1
NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
TATTCAAAAATTGAGAATTTCTGACCACTTAACAAACCCACAGAAAATCCACCCGAGTG
CACTGAGCACGCCAGAAATCAGGTGGCCTCAAAGAGCTGCTCCCACCTGAAGGAGACGCG
CTGCTGCTGCTGTCGTCTGCCTGGCGCCTTGGCCTACAGGGGCCGCGGTTGAGGGTGGG
AGTGGGGGTGCACTGGCCAGCACCTCAGGAGCTGGGGGTGGTGGTGGGGCGGTGGGGGT
GGTGTTAGTACCCCATCTTGAGGTCTGAAACACAAAGTGTGGGGTGTCTAGGGAAGAAG
>2
NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
AAAAGCATTTATGCTACAAATTACTATGGTAATTATGCTACAAATTTATGGTACCATAAA
TTACCATAGTAATTTGTAGCATAAATTTGTACTATGGTACAAATTACATGGGAGAGTGAA
GGTGGGTTAAACATTCAATTAAGAAGTCCACTCAGATTGCAAGAAAAGAGAGAGGA
ATGGAGATGGTAGCACAAAGTCCCTACAATAAAAGTAGATGTTTGTAGATCAGTTCATTT
```

FASTQ

Read 1

@ERR007731.739 IL16_2979:6:1:9:1684/1 ← Read name

CTTGACGACTTGAAAAATGACGAAATCACTAAAAACGTGAAAAATGAGAAATG... ← Sequence

+

BBBCBBBBBBBABBABBBBBBBBABBABBBBBBBBABBABAAAAABBBBBB=>BB... ← Base qualities

Read 2

@ERR007731.740 IL16_2979:6:1:9:1419/1

AAAAAAAAAGATGTCATCAGCACATCAGAAAAGAAGGCAACTTTAAACTTTTC...

+

BBABB/ABABAABABABBABBBAAA>@B@BBAA@4AAA>.>BAA@779:AAA@A...

- | Simple format for raw unaligned sequencing reads
- | Paired-end sequencing: two FASTQ files or one interleaved file

FASTQ

Read 1

```
@ERR007731.739 IL16_2979:6:1:9:1684/1 ← Read name
CTTGACGACTTGAAAAATGACGAAATCACTAAAAACGTGAAAAATGAGAAATG... ← Sequence
+
BBBCBCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC=>BB... ← Base qualities
```

Read 2

```
@ERR007731.740 IL16_2979:6:1:9:1419/1
AAAAAAAAAGATGTCATCAGCACATCAGAAAAGAAGGCAACTTTAAACTTTTC...
+
BBABB/ABABAABABABBABBAAA>@B@BBAA@4AAA>.>BAA@779:AAA@A...
```

- | Simple format for raw unaligned sequencing reads
- | Paired-end sequencing: two FASTQ files or one interleaved file
- | Quality encoded in ASCII characters with decimal codes 33-126
 - | ASCII code of "A" is 65, the corresponding quality is $Q = 65 - 33 = 32$

Base quality encoded as character



! " # \$ % & ' () * + , - . / 0 1 2 3 4 5 6 7 8 9 : ; < = > ? @ A B C D E F G H I J																											
Numeric ASCII value																											
33													47														
Base quality value																											
0													14	(65-33 = 32)													
													32														

Quality = Phred-scaled probability of an error


Quality	Probability of error	Accuracy
10 (Q10)	1 in 10	90%
20 (Q20)	1 in 100	99%
30 (Q30)	1 in 1000	99.9%
40 (Q40)	1 in 10000	99.99%

$$Q = -10 \log_{10} P \quad \dots \quad P = 10^{-Q/10}$$

SAM / BAM: Sequence Alignment/Map format



SAM / BAM: Sequence Alignment/Map format




Flag

Hex	Dec	Flag	Description
0x1	1	PAIRED	paired-end (or multiple-segment) sequencing technology
0x2	2	PROPER_PAIR	each segment properly aligned according to the aligner
0x4	4	UNMAP	segment unmapped
0x8	8	MUNMAP	next segment in the template unmapped
0x10	16	REVERSE	SEQ is reverse complemented
0x20	32	MREVERSE	SEQ of the next segment in the template is reversed
0x40	64	READ1	the first segment in the template
0x80	128	READ2	the last segment in the template
0x100	256	SECONDARY	secondary alignment
0x200	512	QCFAIL	not passing quality controls
0x400	1024	DUP	PCR or optical duplicate
0x800	2048	SUPPLEMENTARY	supplementary alignment

Bit operations made easy

- ```
- samtools flags
 0xa3 163 PAIRED,PROPER_PAIR,MREVERSE,READ2
- python
 0x1 | 0x2 | 0x20 | 0x80 .. 163
 bin(163) .. 10100011
```

## SAM / BAM: Sequence Alignment/Map format



### Insert size

length of the DNA fragment sequenced from both ends by paired-end sequencing:




# CRAM: Reference based Compression

BAM files are too large

- ~1.5-2 bytes per base pair

Increases in disk capacity are being far outstripped by sequencing technologies



Zachary D. Stephens, *et al*, Big Data: Astronomical or Genomical? DOI: 10.1371/journal.pbio.1002195

# CRAM: Reference based Compression

BAM files are too large

- | ~1.5-2 bytes per base pair

Increases in disk capacity are being far outstripped by sequencing technologies

BAM stores all of the data

- | Every read base
- | Every base quality
- | Using a single conventional compression technique for all types of data

```
Reference sequence: ACGTACGTACGTACGTACGTACGTACGTACGTAC
read 1: ACGTACGTACGTACGTACGTGC
read 2: TACGTACGCACGTACGTGCGTA
read 3: CGTACGCACGTACGTACGTACG
read 4: TACGTACGTACGTGCGTACGTA
read 5: CGCACGTACGTACGTACGTACG
read 6: TACGTGCGTACGTACGTAC
```

# CRAM: Reference based Compression

BAM files are too large

- | ~1.5-2 bytes per base pair

Increases in disk capacity are being far outstripped by sequencing technologies

BAM stores all of the data

- | Every read base
- | Every base quality
- | Using a single conventional compression technique for all types of data

```
Reference sequence: ACGTACGTACGTACGTACGTACGTACGTAC
read 1: G.
read 2: C.....
read 3: C.....
read 4: G.....
read 5: ..C.....
read 6: G.....
```

CRAM: in lossless mode 60% of BAM size

- | Reference based compression
- | Controlled loss of quality information
- | Different compression methods for different type of data

# VCF: Variant CallFormat

**VCF header**

```
##fileformat=VCFv4.0
##fileDate=20100707
##source=VCFtools
##reference=NCBI36
##INFO=<ID=AA,Number=1,Type=String,Description="Ancestral Allele">
##INFO=<ID=H2,Number=0,Type=Flag,Description="HapMap2 membership">
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=GQ,Number=1,Type=Integer,Description="Genotype Quality (phred score)">
##FORMAT=<ID=GL,Number=3,Type=Float,Description="Likelihoods for RR,RA,AA genotypes (R=ref,A=alt)">
##FORMAT=<ID=DP,Number=1,Type=Integer,Description="Read Depth">
##ALT=<ID=DEL,Description="Deletion">
##INFO=<ID=SVTYPE,Number=1,Type=String,Description="Type of structural variant">
##INFO=<ID=END,Number=1,Type=Integer,Description="End position of the variant">
```

**Body**

| #CHROM | POS | ID  | REF | ALT   | QUAL | FILTER | INFO               | FORMAT   | SAMPLE1  | SAMPLE2 |
|--------|-----|-----|-----|-------|------|--------|--------------------|----------|----------|---------|
| 1      | 1   | .   | ACG | A,AT  | .    | PASS   | .                  | GT:DP    | 1/2:13   | 0/0:29  |
| 1      | 2   | rs1 | C   | T,CT  | .    | PASS   | H2;AA=T            | GT:GQ    | 0 1:100  | 2/2:70  |
| 1      | 5   | .   | A   | G     | .    | PASS   | .                  | GT:GQ    | 1 0:77   | 1/1:95  |
| 1      | 100 | .   | T   | <DEL> | .    | PASS   | SVTYPE=DEL;END=300 | GT:GQ:DP | 1/1:12:3 | 0/0:20  |

**Annotations:**

- Mandatory header lines** (red arrow pointing to ##fileformat=VCFv4.0)
- Optional header lines** (meta-data about the annotations in the VCF body) (blue arrow pointing to ##INFO=...)
- Reference alleles** (GT=0) (blue arrow pointing to 0/0:29)
- Alternate alleles** (GT>0 is an index to the ALT column) (blue arrow pointing to 1/1:95)
- Phased data** (G and C above are on the same chromosome) (blue arrow pointing to 0|1:100)
- Variant types** (blue arrows pointing to ALT column):
  - Deletion (pointing to <DEL>)
  - SNP (pointing to A,AT)
  - Large SV (pointing to T,CT)
  - Insertion (pointing to G)
  - Other event (pointing to SVTYPE=DEL)

File format for storing variation data

- tab-delimited text, parsable by standard UNIX commands
- flexible and user-extensible
- compressed with BGZF (bgzip), indexed with TBI or CSI (tabix)

# VCF anatomy

```
...
##INFO=<ID=DP,Number=1,Type=Integer,Description="Raw read depth">
##INFO=<ID=AF,Number=A,Type=Float,Description="Allele frequency in population">
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=AD,Number=R,Type=Integer,Description="Allelic depths (ref,alt,...)">
...
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT SAMPLE1 SAMPLE2 SAMPLE3
11 24535 . G A 243 PASS DP=221;AF=0.5 GT:AD 0/1:73,15 0/0:48,0 0/1:71,14
```

Row-oriented, tab-delimited file with eight mandatory columns (CHROM-INFO)

# VCF anatomy

```
...
##INFO=<ID=DP,Number=1,Type=Integer,Description="Raw read depth">
##INFO=<ID=AF,Number=A,Type=Float,Description="Allele frequency in population">
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=AD,Number=R,Type=Integer,Description="Allelic depths (ref,alt,...)">
```

| #CHROM | POS   | ID | REF | ALT | QUAL | FILTER | INFO          | FORMAT | SAMPLE1   | SAMPLE2  | SAMPLE3   |
|--------|-------|----|-----|-----|------|--------|---------------|--------|-----------|----------|-----------|
| 11     | 24535 | .  | G   | A   | 243  | PASS   | DP=221;AF=0.5 | GT:AD  | 0/1:73,15 | 0/0:48,0 | 0/1:71,14 |

Genomic coordinates



# VCF anatomy

```
...
##INFO=<ID=DP,Number=1,Type=Integer,Description="Raw read depth">
##INFO=<ID=AF,Number=A,Type=Float,Description="Allele frequency in population">
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=AD,Number=R,Type=Integer,Description="Allelic depths (ref,alt,...)">
...


#CHROM	POS	ID	REF	ALT	QUAL	FILTER	INFO	FORMAT	SAMPLE1	SAMPLE2	SAMPLE3
11	24535	.	G	A	243	PASS	DP=221;AF=0.5	GT:AD	0/1:73,15	0/0:48,0	0/1:71,14


```

Arbitrary string, typically a dbSNP RefSNP id. Dot for missing value.


# VCF anatomy

```
...
##INFO=<ID=DP,Number=1,Type=Integer,Description="Raw read depth">
##INFO=<ID=AF,Number=A,Type=Float,Description="Allele frequency in population">
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=AD,Number=R,Type=Integer,Description="Allelic depths (ref,alt,...)">
...
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT SAMPLE1 SAMPLE2 SAMPLE3
11 24535 . G A 243 PASS DP=221;AF=0.5 GT:AD 0/1:73,15 0/0:48,0 0/1:71,14
```



# VCF anatomy

```
...
##INFO=<ID=DP,Number=1,Type=Integer,Description="Raw read depth">
##INFO=<ID=AF,Number=A,Type=Float,Description="Allele frequency in population">
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=AD,Number=R,Type=Integer,Description="Allelic depths (ref,alt,...)">
...
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT SAMPLE1 SAMPLE2 SAMPLE3
11 24535 . G A 243 PASS DP=221;AF=0.5 GT:AD 0/1:73,15 0/0:48,0 0/1:71,14
```



# VCF anatomy

```
...
##INFO=<ID=DP,Number=1,Type=Integer,Description="Raw read depth">
##INFO=<ID=AF,Number=A,Type=Float,Description="Allele frequency in population">
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=AD,Number=R,Type=Integer,Description="Allelic depths (ref,alt,...)">
...
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT SAMPLE1 SAMPLE2 SAMPLE3
11 24535 . G A 243 PASS DP=221;AF=0.5 GT:AD 0/1:73,15 0/0:48,0 0/1:71,14
```

Although in theory phred-scaled probability, don't expect truly probabilistic interpretation in practice.

# VCF anatomy

```
...
##INFO=<ID=DP,Number=1,Type=Integer,Description="Raw read depth">
##INFO=<ID=AF,Number=A,Type=Float,Description="Allele frequency in population">
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=AD,Number=R,Type=Integer,Description="Allelic depths (ref,alt,...)">
...
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT SAMPLE1 SAMPLE2 SAMPLE3
11 24535 . G A 243 PASS DP=221;AF=0.5 GT:AD 0/1:73,15 0/0:48,0 0/1:71,14
```

Soft-filter variants with e.g. low quality, low depth, etc.

# VCF anatomy

```
...
##INFO=<ID=DP,Number=1,Type=Integer,Description="Raw read depth">
##INFO=<ID=AF,Number=A,Type=Float,Description="Allele frequency in population">
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=AD,Number=R,Type=Integer,Description="Allelic depths (ref,alt,...)">
...
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT SAMPLE1 SAMPLE2 SAMPLE3
11 24535 . G A 243 PASS DP=221;AF=0.5 GT:AD 0/1:73,15 0/0:48,0 0/1:71,14
```

Per-site annotations. Here **DP** is the cumulative read depth across all samples and **AF** allele frequency of the allele in general population.

# VCF vs BCF

VCFs can be very big

- compressed VCF with 3781 samples, human data:
  - 54 GB for chromosome 1
  - 680 GB whole genome

VCFs can be slow to parse

- text conversion is slow
- main bottleneck: FORMAT fields


```
##fileformat=VCFv4.0
##fileDate=20100707
##source=VCFtools
##ALT=<ID=DEL,Description="Deletion">
##INFO=<ID=END,Number=1,Type=Integer,Description="End position of the variant">
##CHROM POS ID REF ALT QUAL FILTER INFO FORMAT SAMPLE1 SAMPLE2
1 3 . A G . PASS AC=67;AN=5400;DP=2809 GT:PL:DP:GQ 1/1:0,9,73:26:22 0/0:0,9,73:13:31 0/0:0,9,73:48:99 1/0:255,0,75:32:15 1/0:255,0,75:32:15
1 4 . A T . PASS AC=15;AN=6800;DP=6056 GT:PL:DP:GQ 0/0:0,9,73:13:31 1/0:255,0,75:32:15 0/0:0,2,80:14:90 1/1:0,9,73:26:22 0/0:0,9,73:13:31
1 5 . C T . PASS AC=20;AN=6701;DP=5234 GT:PL:DP:GQ 1/0:255,0,75:32:15 0/0:0,2,170:14:90 1/1:0,9,73:13:31 0/0:0,6,50:13:80 0/0:0,2,80:14:90
1 6 . A G . PASS AC=67;AN=5400;DP=2809 GT:PL:DP:GQ 1/1:0,9,73:26:22 0/0:0,9,73:13:31 0/0:0,9,73:48:99 1/0:255,0,75:32:15 1/0:255,0,75:32:15
1 7 . A T . PASS AC=15;AN=6800;DP=6056 GT:PL:DP:GQ 0/0:0,9,73:13:31 1/0:255,0,75:32:15 0/0:0,2,80:14:90 1/1:0,9,73:26:22 0/0:0,9,73:13:31
```

BCF

- binary representation of VCF
- fields rearranged for fast access

| #CHROM | POS | ID | REF | ALT | QUAL | FILTER | INFO         | FORMAT                 | SAMPLE1                            | SAMPLE2           | SAMPLE3           | SAMPLE4            | SAMPLE5            |
|--------|-----|----|-----|-----|------|--------|--------------|------------------------|------------------------------------|-------------------|-------------------|--------------------|--------------------|
| 1      | 6   | .  | A   | G   | .    | PASS   | AC=67;AN=540 | GT:PL:DP:GQ            | 1/1:0,9,73:26:22                   | 0/0:0,9,73:13:31  | 0/0:0,9,73:48:99  | 1/0:255,0,75:32:15 | 1/0:255,0,75:32:15 |
| 1      | 6   | .  | A   | G   | .    | PASS   | AC=67;AN=540 | GT:1/1:0/0:0/0:1/0:1/0 | PL:0,9,73:0,9,73:255,0,75:255,0,75 | DP:26:13:48:32:32 | GQ:22:31:99:15:15 |                    |                    |

# File formats summary





# File formats summary

